


High Quality Dental Instruments



SURGIFUSION
INTERNATIONAL



✉ info@surgifusionintl.com

☎ +92 313 4260504

🌐 www.surgifusionintl.com

📍 **Address:**
Roras Road new Miana Pura,
Sialkot 51310, PAKISTAN



Solutions for dental surgeons

We are a trusted Manufacturer & Exporter of Dental Instruments, dedicated to delivering high-quality, precision-engineered instruments for dental professionals, clinics, hospitals, and distributors worldwide. With a strong focus on quality, innovation, and customer satisfaction, we manufacture a comprehensive range of dental instruments designed to meet the demanding requirements of modern dentistry.

Our product range includes dental examination instruments, extraction forceps, elevators, orthodontic instruments, periodontal instruments, restorative instruments, surgical instruments, and other specialized dental tools. Every instrument is manufactured using carefully selected materials and advanced production techniques to ensure excellent performance, durability, precision, and long-lasting reliability.

Our skilled workforce, modern manufacturing facilities, and strict quality control procedures enable us to maintain consistent standards throughout every stage of production. From raw material selection and precision manufacturing to finishing, inspection, and packaging, we pay close attention to every detail.

As an experienced exporter, we proudly serve customers and business partners in international markets. We offer reliable service, competitive pricing, flexible manufacturing solutions, and timely delivery to meet the specific requirements of our global clients.

Our mission is to build long-term partnerships by providing dependable dental instruments that professionals can trust. We continuously strive to improve our products and services while maintaining our commitment to quality, precision, and customer satisfaction.





Care for Dental Instruments

“The cost of maintenance is lesser than the cost of new instruments.”

The Quality, To Begin With

Every Surgifusion International passes through a system of strict quality control after each manufacturing step, guaranteeing the highest quality only. The quality control is both the inherent one and the one complying with international standards.

The Steel, To Be Sure

We procure our own special steel with a higher percentage of carbon and chromium. Instruments fabricated out of this steel guarantee a longer life and durable functionality, ensuring satisfaction to users.

Experienced Craftsmen, To Create Master Pieces

All instruments are crafted by hand. Be it forging, machining, grinding and sharpening, all these processes are performed by our trained and experienced craftsmen. During manufacturing, every single instrument passes through numerous steps and every step is followed by an inspection process. We do understand what depends on the high quality of our instruments and to this end we allocate our best resources.

Intended And Specified Use, The Long Life

Always use the instrument for the function it has been designed for. Dental instruments are ultra fine precision instruments. A misuse can result in severe damage to the working points. Increase the life of instrument by its proper and correct use only. A pair of scissors is meant for cutting soft tissues only. It should not be used to lift bones or for exploring and scraping purpose.

Disinfection, Cleaning And Sterilization

Always clean the instruments before sterilization. If the instruments are sterilized without cleaning, the stains on them can be burnt into the instruments leaving permanent marks. This is because of very high temperature during sterilization.

Before cleaning, instruments should be disinfected using a suitable disinfectant preferably the one not containing Phenol.

Manual cleaning should begin soon after the surgeries are over. Do not delay cleaning. Use a soft brush and not a steel wool or drill brush to clean the instruments to remove stubborn impurities. Instruments should always be cleaned in a non-corrosive, natural cleaning agent with little foaming. The instruments should then be rinsed in distilled water.

The use of ultrasonic cleaners and thermo-disinfectant is also becoming common. Here again, care must be taken in the selection of cleaning agents, detergents and disinfectants. Please ensure that they do not damage the metal. Always remember these methods can clean but do not sterilize!



Some Important Care Tips:

- Before sterilization, always inspect your equipment for any dirt remains or any organic or mineral deposits. The extreme heat can cause these deposits to be permanently burnt into the instrument, increasing chances of corrosion.
- Never mix stainless steel with other metals like aluminum, brass or copper. This can cause electrochemical changes and also produce etching or staining.
- Use distilled or demineralized water only.
- Always follow the manufacturer's instructions precisely when it comes to the temperatures and times involved for autoclaving, dry heating or germicidal solvents.
- Dry each instrument separately and completely. In case of remaining moisture, dry using a sterile cloth.
- Always use physiological oil and other suitable lubricants to keep instruments well oiled. This is especially required for instruments comprising many parts like scissors, pliers and gouges. A regular use of oil prevents stiffness and corrosion and results in smooth functioning.
- Hinged instruments like forceps and needle holders should always be opened for cleaning and sterilization.
- Etching or engraving instruments for personal identification can damage the instrument itself and cause damage to other instruments as well. Damage to the polished surface can always lead to corrosion.
- Although stainless steel has outstanding corrosion resistant properties, contact with certain chemicals must be avoided.





SURGIFUSION
INTERNATIONAL

Table of Contents

Chapter 1	01	Chapter 5	25
Dental Extraction Forceps Zahnzangen Forceps de extracción Daviers Pinze per estrazioni		Root Elevators Wurzelherber Elevadores para raíces Elévateurs à racines Leve per radici	
English Pattern Englische Form Forma inglesa Forme anglaise Forma inglese		Chapter 6	33
Chapter 2	11	Amalgam Instruments Amalgaminstrumente Instrumentos para amalgama Instruments à amalgame Strumenti per amalgma	
Dental Extraction Forceps Zahnzangen Forceps de extracción Daviers Pinze per estrazioni		Chapter 7	45
With Alpha Handle Mit Alfa-Griff Con mango Alfa Manche Alfa		Scalers, Curettes Zahnreiniger, Küretten Limpiadores de sarro, Curetas Instruments à détartrer, Curettes Scaler, Curette	
Chapter 3	12	Chapter 8	49
Dental Extraction Forceps Zahnzangen Forceps de extracción Daviers Pinze per estrazioni		Periodontia Instruments Instrumente Zur paradontologie Instrumentos para periodoncia Instruments pour la parodontologie Strumenti per parodontologia	
For children Kinderzahnzangen Para niños Pour enfants Per bambini		Chapter 9	57
Chapter 4	16	Excavators, Probes & Explorers, Nerve instruments Exkavatoren, Zahnsonden, Nervinstrumente Excavadores, sondas y exploradores Instrumentos para nervios Excavateurs, Sondes, et explorateurs Escavatori, Sonde, Strumenti per nervi	
Dental Extraction Forceps Zahnzangen Forceps de extracción Daviers Pinze per estrazioni		Chapter 10	67
American Pattern Amerikanische form Forma americana Forme américaine Forma americana		Wax and Crown Instruments, Matrix Instruments Wachs-kroneninstrumente, Matrizen-instrumente Instrumentos para cera, coronas, instrumentos para matrices Instruments pour couronnes et matrices Strumenti per cera e per corone	

Chapter 11 75

Impression Trays

Abdrucklöffel
Portamoldes
Porte-empreintes
Porta-impronte

Chapter 12 87

Pliers For Orthodontics

Zangen für Kieferorthopädie
Pinzas para orthodoncia
Pinces pour orthodontie
Pinze per Ortodonzia

Chapter 13 97

Knives, Scalpels

Messer, Skalpel
Bisturíes
Bistouris
Bisturi

Chapter 14 99

Scissors

Scheren
Tijeras
Ciseaux
Forbici

Chapter 15 102

Forceps

Pinzetten
Pinzas
Pinces
Pinzette

Hemostatic Forceps

Arterienklemmen
Pinzas hemostáticas
Pinces hémostatiques
Pinze emostatiche

Chapter 16 107

Needle holders

Nadel halter
Portaagujas
Porte-aiguilles
Porta aghi

Chapter 17 109

Retractor

Wundhaken
Separadores
Ecarteurs
Retrattori

Chapter 18 112

Bone Cutting Forceps, Bone Rongeurs

Hohlemeißelzangen,
Knochen Splitterzangen
Pinzas Cortantes para huesos
Pinces à os coupantes,
Pinces-gouges
Scalpellini e pinze alveolotome

Tongue Depressors, Mouth Gags

Zungenspatel, Mundsperrer
Bajalenguas, Abrebocas
Abaisse-langues, Ouvre-bouches
Spatola per lingu, a Abribocca

Bone Chisels, Curettes

Flachmeißel, Scharfe Löffel
Cinceles, Cucharillas Cortantes
Ciseaux-burins, Curettes tranchantes
ScalPELLi, Cucchiaini per osso

Chapter 19 119

Sterilization, Asepsis

Sterilisation, Asepsis
Esterilización, Asepsia
Stérilisation, Asepsis
Sterilizzazione, Asepsi

Chapter 20 125

Surgery

Chirurgie
Chirguia
Chirurgie
Chirurgia

Chapter 21 131

Dental instruments sets

Satz Instrumente
Sistema de instrumentos
Ensemble d' instruments
Insieme degli strumenti

CH 01-04

- CH 01** English Pattern
Englische Form
Forma inglesa
Forma anglaise
Forma inglese
- CH 02** With Alpha Handle
Mit Alfa-Griff
Con mango Alfa
Manche Alfa
Maniglia dell'alfa
- CH 03** For children
Kinderzahnzangen
Para niños
Pour enfants
Per bambini
- CH 04** American Pattern
Amerikanische form
Forma americana
Forme américaine
Forma americana

Dental Extraction Forceps
Forceps de extracción
Zahnzangen
Daviers
Pinze per estrazioni





SURGIFUSION
INTERNATIONAL

Extracting Forceps, English Pattern

Zahnzangen, englische form

Forceps de extracción, forma inglesa

Daviers, forma anglaise

Pinze per estrazioni, forma inglese

*All Extracting Dental Forceps are also available in **Non-Slip Grip Jaws**.
Please see page 11 for ordering details.*



SI- 01-0101
Fig.1



SI- 01-0102
Fig.2



SI- 01-0103
Fig.3



SI- 01-0104
Fig.4



SI- 01-0107
Fig.7



SI- 01-0108
Fig.8



SI- 01-0113
Fig.13



SI- 01-0117
Fig. 17



SI- 01-0118
Fig. 18



SI- 01-0218
Fig. 18 A



SI- 01-0219
Fig. 19



SI- 01-0220
Fig. 20



SI- 01-0221
Fig. 21



SI- 01-0222
Fig. 22

*All Extracting Dental Forceps are also available in **Non-Slip Grip Jaws**.
Please see page 11 for ordering details.*




SI- 01-0223
Fig. 23



SI- 01-0224
Fig. 24



SI- 01-0229
Fig. 29



SI- 01-0329
Fig. 29 N



SI- 01-0330
Fig. 30



SI- 01-0331
Fig. 31



SI- 01-0333
Fig. 33



SI- 01-0433
Fig. 33 A

Extracting Forceps, English Pattern

Zahnzangen, englische form

Forceps de extracción, forma inglesa

Daviers, forma anglaise


Pinze per estrazioni, forma inglese



SI- 01-0533
Fig. 33 C



SI- 01-0633
Fig. 33 M




SI- 01-0641
Fig. 41



SI- 01-0644
Fig. 44




SI- 01-0645
Fig. 45




SI- 01-0651
Fig. 51


*All Extracting Dental Forceps are also available in **Non-Slip Grip Jaws**.
Please see page 11 for ordering details.*




SI- 01-0751
Fig. 51 A




SI- 01-0851
Fig. 51 C




Mod. Pro
SI- 01-0951
Fig. 51 M



SI- 01-0952
Fig. 52



SI- 01-0955
Fig. 55



SI- 01-0956
Fig. 56


Extracting Forceps, English Pattern

Zahnzangen, englische form


Forceps de extracción, forma inglesa

Daviers, forma anglaise

Pinze per estrazioni, forma inglese




SI- 01-0967
Fig. 67



SI- 01-1067
Fig. 67 A



SI- 01-1073
Fig. 73



SI- 01-1074
Fig. 74



SI- 01-1174
Fig. 74 M




SI- 01-1274
Fig. 74 N



SI- 01-1275
Fig. 75

All Extracting Dental Forceps are also available in **Non-Slip Grip Jaws**.
 Please see page 11 for ordering details.



1/1
Lowrence Read
SI- 01-1276
 Fig. 76



Lowrence Read
SI- 01-1376
 Fig. 76 N



SI- 01-1379
 Fig. 79



SI- 01-1479
 Fig. 79 A




1/1
SI- 01-1486
 Fig. 86 A



SI- 01-1586
 Fig. 86 B



SI- 01-1686
 Fig. 86 C



SI- 01-1687
 Fig. 87


Extracting Forceps, English Pattern

Zahnzangen, englische form

Forceps de extraacción, forma inglesa


Daviers, forma anglaise

Pinze per estrazioni, forma inglese




SI- 01-1689

Fig. 89




SI- 01-1690

Fig. 90



SI- 01-1694

Fig. 94



SI- 01-1695


Fig. 95



SI- 01-1699


Fig. 99

*All Extracting Dental Forceps are also available in **Non-Slip Grip Jaws**.
Please see page next for ordering details.*




1/1


Mead
SI- 01-1701
Fig. MD 1




Mead
SI- 01-1702
Fig. MD 2



Mead
SI- 01-1703
Fig. MD 3



Mead
SI- 01-1704
Fig. MD 4



Alpha handle

All preceding extracting forceps in English pattern are also available with **Alfa handle**.



Just replace first two article numbers, “**01-**” with “**02-**”. Please refer to column “C” below.

Non-Slip Grip Jaws



All preceding extracting forceps in English pattern are also available in **Non-Slip Grip Jaws**.

Please add “NS” with existing article numbers. e.g **SI- 01-0101NS**. Please refer to column B.





SI- 03-1839
Fig. 39
For upper molars, either side



SI- 03-1939
Fig. 39 A
For upper molars, either side



SI- 03-2039
Fig. 39 L
For left upper molar




03-1810
Instrument set of 10




SI- 03-2139
Fig. 39 R
For right upper molar



SI- 03-2140
Fig. 40
For lower molars, either side




SI- 03-2151
Fig. 51 S
For upper roots, either side




Klein
SI- 03-2203
Fig. 3
For upper molars



Klein
SI- 03-2205
Fig. 5
For lower incisors



Klein
SI- 03-2206
Fig. 6
For lower molars



Klein
SI- 03-2207
Fig. 7
For lower roots




Klein
SI- 03-3137
Fig. 137
For upper incisors & canines



Klein
SI- 03-3139
Fig. 139
For upper bicusps




Klein
SI- 03-3151
Fig. 51 S
For upper roots





Standard
SI- 04-0101
Fig. 1
For upper centrals,
laterals and cuspids




Henahan
SI- 04-0201
Fig. 1 A
For upper incisors



Physick
SI- 04-0205
Fig. 5
For lower third molars
either side



SI- 04-0206
Fig. 6
Splitting or separating
forceps



Woodward
SI- 04-0203
Fig. 3 FH
For lower molars



SI-04-0310 SI-04-0410
Fig. 10 H / Fig. 10 S
For upper molars, either side



SI- 04-0310
Fig. 10 H




SI- 04-0410
Fig. 10 S



Harris
SI- 04-0415
Fig. 15




SI- 04-0416
Fig. 16
For lower molars







SI- 04-0723
Fig. 23
For lower molars, either side



SI- 04-0724
Fig. 24
For upper molars, either side



SI- 04-0727
Fig. 27



SI- 04-0732 **SI- 04-0832**
Fig. 32 Fig. 32 A
Alveolar forceps for upper molars
and bicuspsids, strong beak

Parmly
SI-04-0732 - SI-04-0832




SI- 04-0853
Fig. 53 L
For upper
molars left side




SI- 04-0953
Fig. 53 R
For upper
molars right side



SI- 04-0962
 Fig. 62
 Universal forceps for upper
 teeth & roots, also for children



SI- 04-0965
 Fig. 65
 For upper incisors
 and roots, long
 slender beak




Tomes
SI- 04-0969
 Fig. 69
 For fragments and
 small roots upper & lower



Hull
SI- 04-1101
Fig .101
Universal forceps for
bicuspid and milk teeth

SI- 04-1103
Fig .103
For lower bicuspid
incisors, roots, either side



Nevius
SI- 04-1188
Fig .88 L
For left upper molar

Nevius
SI- 04-1288
Fig .88 R
For right upper molar



Kells
SI- 04-1299
Fig .99 A
For upper bicuspid,
cuspid and incisors



Kells
SI- 04-1399
Fig .99 C
For upper bicuspid,
cuspid and incisors



Nevius
SI- 04-2110
Fig .110
For upper cuspid and bicuspid



Nevius
SI- 04-3110
Fig .110
For upper cuspid and bicuspid



Cryer
SI- 04-4150
Fig. 150



Cryer
SI- 04-5150
Fig. 150 A



Cryer
SI- 04-6150
Fig. 150 S




Cryer
SI- 04-6151
Fig. 151




Cryer
SI- 04-7151
Fig. 151 A



Cryer
SI- 04-8151
Fig. 151 S




SI- 04-8203
Fig. 203




SI- 04-8210 - 04-9210

SI- 04-8210
Fig. .210 H


SI- 04-9210
Fig. .210 S



SI- 04-9217
Fig. .217



SI- 04-9222
Fig. .222



SI- 04-9286
Fig. .286

CH 05

Root Elevators
Wurzelheber
Elevadores para raíces
Elévateurs à racines
Leve per radici



Choice of Handles for Root Elevators.


If you wish to order a different type (**A, B, C, D, E, F and G**) of handle for a particular product, please mention that type with the code number of that product.




Example : SI-05-0701E

If you wish to order a particular type (**A, B, C, D, E, F and G**) of handle with all instruments, please mention that type in your purchase order.


**Example : Please use handle C for all
Scalers.**









SI-05-0934 - SI-05-2304
13.5cm/ 5 1/4"





Winter

SI-05-2801 - SI-05-3314



Walter F. Barry

SI-05-4320 - SI-05-4371



Fig. 1

Krallengnber
SI- 05-4401




Fig. 2

Krallengnber
SI- 05-4402



Lecluse
SI- 05-4501




1/1

SI- 05-4601

SI- 05-4602

Pott
SI-05-4604 - SI-05-4602




CH 06

Amalgam Instruments
Amalgaminstrumente
Instrumentos para amalgama
Instruments à amalgame
Strumenti per amalgma



**Choice of Handles for Scalars, Amalgam Instruments,
Filling Instruments, Curettes, Probes, Knives, Explorers and Excavators**





If you wish to order a different type (**A, B, C, D, E and F**) of handle for a particular products, please mention that type with the code number of that product.

Example : SI-06-1001C





If you wish to order a particular type (**A, B, C, D, E and F**) of handle with all instruments, please mention that type in your purchase order.

Example : Please use handle C for all Scalers.





SI-06-0101 - SI-06-0315





Amalgam well
SI- 06-0935
3,5 cm Ø
3,5 cm h.





SI- 06-1201
Fig.1

SI- 06-1202
Fig.2

SI- 06-1303
Fig.3

SI- 06-1403
Fig.3 A


SI- 06-1404
Fig.4

SI- 06-1504
Fig.4 A

SI- 06-1544
Fig.4 R

Mod. PHF


SI- 06-1201 -
SI- 06-1544
16cm/6 1/4"



SI- 06-1604 **SI- 06-1606**
Fig.4 Fig.6

SI- 06-1700 **SI- 06-1701** **SI- 06-1702** **SI- 06-1703**
Fig.0 Fig.1 Fig.2 Fig.3

Williams



SI- 06-1800 Fig.0 SI-06-1801 Fig.1 SI- 06-1802 Fig.2 SI- 06-1903 Fig.1 SI- 06-1904 Fig.3 SI- 06-1905 Fig.2 SI- 06-1906 Fig.3

White

Woodson


SI- 06-1800 -
SI- 06-2205




SI- 06-2103 Fig.3 SI- 06-2104 Fig.4 SI- 06-2105 Fig.5 SI- 06-2201 Fig.1 SI- 06-2202 Fig.2 SI- 06-2203 Fig.3 SI- 06-2204 Fig.4 SI- 06-2205 Fig.5

White


Ash




Please see page 34 for different handle styles.




SI-06-2713 - SI-06-3037




SI- 06-2713
Fig.13




SI- 06-2714
Fig.14




SI- 06-2715
Fig.15




SI- 06-2716
Fig.16




SI- 06-2826
Fig.26




SI- 06-2827
Fig.27




SI- 06-2828
Fig.28




SI- 06-2917
Fig.17




SI- 06-2918
Fig.18




SI- 06-2919
Fig.19



SI- 06-3035
Fig.35



SI- 06-3036
Fig.36



SI- 06-3037
Fig.37

Burnishers

Polierinstrumente

Bruñidores

Brunissoirs et Couteaux à amalgame

Strumenti per finitura



1/1
White
SI-06-3101 SI-06-3105



SI- 06-3101



SI- 06-3102




SI- 06-3103



SI- 06-3104



SI- 06-3105




SI- 06-3201
Fig.1



SI- 06-3202
Fig.2



SI- 06-3203
Fig.3



SI- 06-3301




SI- 06-3302




SI- 06-3304

SI-06-3201 - SI-06-3203



1/1

SI-06-3501




SI-06-3502



SI-06-3503




SI-06-3506



SI- 07-0134
Fig.33 ½

SI- 07-0135
Fig.34 ½



Darby-Perry

SI- 07-0201
Fig.1

SI- 07-0204
Fig.4

SI- 07-0205
Fig.5

SI- 07-0206
Fig.6

SI- 07-0211
Fig.11



SI- 07-0301
Fig.1



SI- 07-0302
Fig.2



SI- 07-0303
Fig.3



SI- 07-0322
Fig.22



SI- 07-0333
Fig.33



SI- 07-0344
Fig.44

SI-07-0134 - SI-07-0618
15.5 cm/ 6"



SI- 07-0412
Fig.12



SI- 07-0533
Fig.33



SI- 07-0534
Fig.34



SI- 07-0535
Fig.35

Jaquette




SI- 07-0617
Fig.17



SI- 07-0618
Fig.18

McCall

Please see page 34 for different handle styles.






Younger-Good
SI- 07-0978
Fig. 7/8




Barnhart
SI- 07-1012
Fig.   



Barnhart
SI- 07-1056
Fig. 5/6




Indiana University (Anterior)
SI- 07-1134
Fig. 13/14+




Indiana University (Posterior)
SI- 07-1178
Fig. 17/18+



Bunting
SI- 07-1212
Fig.   +



Bunting
SI- 07-1256
Fig. 5/6



Bates
SI- 07-1334
Fig. 135/145



Rule
SI- 07-1434
Fig.   



SI- 07-1441
Fig. 411/412




SI- 07-1541
Fig. 413/414




Mallery
SI- 07-1612
Fig.   



Goldman
SI- 07-1632
Fig. 31/32



SI- 07-1711
Fig. 1





Anterior
SI- 07-1812
Fig.   

Please see page 34 for different handle styles.

CH 08

Periodontia Instruments
Instrumente Zur paradontologie
Instrumentos para periodoncia
Instruments pour la parodontologie
Strumenti per parodontologia






SI- 08-0100
Fig. 8



SI- 08-0101
Fig. W



SI- 08-0108
Fig. 11



SI- 08-0110
Fig. 0



SI- 08-0111
Fig. 10



SI- 08-0112
Fig. WHO




SI- 08-0201




SI- 08-0301




SI- 08-0302







Please see page 34 for different handle styles.



Goldman-Fox Nabers Cross Prichard
SI- 08-1507 SI- 08-1508 SI- 08-1509 SI- 08-1510 SI- 08-1511 SI- 08-1521 SI- 08-1600 SI- 08-1601 SI- 08-1700



Please see page 34 for different handle styles.



Ochsenbein
SI- 08-2201
Fig. 1

Ochsenbein
SI- 08-2202
Fig. 2

Ochsenbein
SI- 08-2303
Fig. 3

Ochsenbein
SI- 08-2304
Fig. 4

Fedi
SI- 08-2501
Fig. 1

Fedi
SI- 08-2502
Fig. 2

Fedi
SI- 08-2503
Fig. 3



SI- 08-2612
Fig. 1/2

SI- 08-2634
Fig. 3/4

SI- 08-2656
Fig. 5/6

Gracey
SI- 08-2678
Fig. 7/8

SI- 08-2690
Fig. 9/10

SI- 08-2692
Fig. 11/12


SI- 08-2694
Fig. 13/14

Please see page 34 for different handle styles.

CH 09


**Excavators, Probes &
Explorers, Nerve instruments**
Exkavatoren, Zahnsonden,
Nervinstrumente
Excavadores, sondas y
exploradores, Instrumentos
para nervios
Excavateurs, Sondes, et
explorateurs
Escavatori, Sonde, Strumenti
per nervi






SI- 09-0105 SI- 09-0113 SI- 09-0114 SI- 09-0117 SI- 09-0118 SI- 09-0131 SI- 09-0133 SI- 09-0135

Please see page 34 for different handle styles.




SI- 09-0212 -
SI- 09-0279

SI- 09-0305 -
SI- 09-0320




SI- 09-0212 Fig. 12/13 SI- 09-0214 Fig. 14/15 SI- 09-0217 Fig. 171/172 SI- 09-0273 Fig. 173/174 SI- 09-0275 Fig. 175/176 SI- 09-0277 Fig. 177/178 SI- 09-0279 Fig. 179/180



SI- 09-0305 Fig. 5 SI- 09-0314 Fig. 14 SI- 09-0317 Fig. 17 SI- 09-0318 Fig. 18 SI- 09-0319 Fig. 19 SI- 09-0320 Fig. 20

White





**SI- 09-0503-
SI- 09-0525**



SI- 09-0503
Fig. 3/4



SI- 09-0505
Fig. 5/6



SI- 09-0507
Fig. 7/8



SI- 09-0509
Fig. 9/10



SI- 09-0511
Fig. 11/12




SI- 09-0513
Fig. 13/14



SI- 09-0515
Fig. 15/16



SI- 09-0517
Fig 17/18



SI- 09-0519
Fig 19/20





SI- 09-0521
Fig 21/22



SI- 09-0523
Fig 23/24



SI- 09-0525
Fig 25/26



Please see page 34 for different handle styles.



I SI- 09-0803	SI- 09-0806	SI- 09-0808	SI- 09-0809	SI- 09-0810
V SI- 09-0903	SI- 09-0906	SI- 09-0908	SI- 09-0909	SI- 09-0910
Fig. 3	Fig. 6	Fig. 8	Fig. 8 a	Fig. 9



I SI- 09-0817	SI- 09-0823	SI- 09-0833	SI- 09-0835
V SI- 09-0917	SI- 09-0923	SI- 09-0933	SI- 09-0935
Fig. 17	Fig. 23	Fig. 33	Fig. 35



I SI- 09-0842	SI- 09-0843	SI- 09-0844
V SI- 09-0942	SI- 09-0943	SI- 09-0944
Fig. 42	Fig. 43	Fig. 44



I SI- 09-0854	SI- 09-0856	SI- 09-0861
V SI- 09-0954	SI- 09-0956	SI- 09-0961
Fig. 54	Fig. 56	Fig. 61



I V
Grip Types


Nerve-canal and root-canal condensers

Nervkanal-und wurzelkanal-stopfer

Tapón para el canal del nervio y el canal de la raíz

Obturbateur de canal nerveux et de canal radiculaire

Condensatori per nervi e canali radicolari



1/1
SI- 09-1025
0,25 mm


SI- 09-1040
0,40 mm

SI- 09-1060
0,60 mm

SI- 09-1073
Fig.3

SI- 09-1080
Fig.0

SI- 09-1090
Fig.00



Files / Black Cutting Instruments

Feilen/ Schneidende instrumente nach Black

Limas para raices/ instrumentos cortantes de Black

Limes/ instruments tranchantes de Black


Lime/ Strumenti da taglio di Black




Hirschfeld
SI- 09-1237
Fig. 3/7



Hirschfeld
SI- 09-1211
Fig. 5/11



Black
SI-09-1351 - SI-09-1379



SI- 09-1351


Fig. 51 / Fig. 52
a) 15 15
b) 8 8
c) 12 L 12 R

SI- 09-1353

Fig. 53 / Fig. 54
a) 10 10
b) 6 6
c) 12 L 12 R

SI- 09-1363

Fig. 63 / Fig. 64
a) 15 15
b) 8 8
c) 12 L 12 R



SI- 09-1365

Fig. 65 / Fig. 66
a) 10 10
b) 8 6
c) 12 L 12 R

SI- 09-1377

Fig. 77 / Fig. 78
a) 15 15
b) 8 8
c) 12 L 12 R

SI- 09-1379

Fig. 79 / Fig. 80
a) 15 15
b) 8 8
c) 12 L 12 R

a) Width of blade in 1/10 mm
b) Length of blade from the centre of angle in mm
c) Angle of blade to its handle in centigrade

CH 10


**Wax and Crown Instruments,
Matrix Instruments**


**Wachs- kroneninstrumente,
Matrizen-instrumente**

**Instrumentos para cera, coronas,
instrumentos para matrices
Instruments pour couronnes
et matrices**

**Strumenti per cera e per corone,
Strumenti per matrici**







Wax, Porcelain and Modelling Carvers

Wachs-porzellan und modellierinstrumente

Talladores y espátulas para cera, porcelana y para modelar

Instruments à cire et à modeler

Strumenti per cera e modellazione



Le Cron
SI- 10-0601



Le Cron
SI- 10-0602



Fehnenstock
SI- 10-0713 - 13 cm/ 5"
SI- 10-0717 - 17 cm/ 6 ¾"



Ermert
SI- 10-0801



Beale
SI- 10-0807



Gisy-Wild
SI- 10-0815



Miller
SI- 10-0915



Sachs
SI- 10-1001



Mershon
SI- 10-1101



SI- 10-1201



1/1

SI- 10-1301



SI- 10-1302




SI- 10-1303



SI- 10-1304



SI- 10-1305




Treyemann
SI- 10-1306



Schroeder
SI-10-1300
with interchangeable
points 10-1301-10-1305



SI- 10-1307



Bohm
SI- 10-1511
11.5 cm/ 4½"

Bohm
SI- 10-1514
11.5 cm/ 4½"



SI- 10-1401 SI- 10-1402 SI- 10-1403



SI-10-1401 - SI-10-1403



SI- 10-1600 SI-10-1601 SI-10-1602
(Handle for
DS-10-1600-RS-10-1602



SI-10-1800

For Molars

For Molars and Premolars



1/1


Tofflemire-Universal
SI- 10-1905




SI- 10-1907
Fig. 7



SI- 10-1908
Fig. 8




SI- 10-1909
Fig. 9



SI- 10-1910
Fig. 10



SI- 10-1911
Fig. 1



SI- 10-1912
Fig. 12




SI- 10-1913
Fig. 13

(Package of 12 bands)




1/1


Ivory
SI- 10-2005 - 5 mm
SI- 10-2006 - 6 mm
SI- 10-2007 - 7 mm




Tofflemire
SI- 10-2101
Fig. 1



Tofflemire
SI- 10-2102
Fig. 2



Tofflemire
SI- 10-2103
Fig. 3




Tofflemire
SI- 10-2104
Fig. 4




1/1
Ivory
SI- 10-2201
Fig. 1
(Without shield)

Ivory
SI- 10-2212
Fig. 12
(With shield)

Ivory
SI- 10-2312
Fig. 12 a



SI-10-2401-SI-10-2403



SI-10-2404-SI-10-2406



Ainsworth
SI- 10-2516
16 cm/ 6 1/4"



SI- 10-2515
15cm/6"



Stockes
SI- 10-2617
17 cm/ 6 3/4"



Brewer
SI- 10-2717
17 cm/ 6 3/4"



SI- 10-2716
16 cm/ 6 1/4"

CH 11


Impression Trays
Abdrucklöffler
Portamoldes
Porte-empreintes
Porta-impronte



Stainless Steel



Stainless steel “Solid” impression trays with retention rim.



set
10

Set of 10 pieces.

SI-11-2003

For Upper Jaw






For Lower Jaw



		mm.	mm.
SI-11-2051- U1	XL	81	63
SI-11-2052- U2	L	74	60
SI-11-2053- U3	M	71	55
SI-11-2054- U4	NM	67	53
SI-11-2055- U5	S	65	51

		mm.	mm.
SI-11-2061- L1	XL	80	63
SI-11-2062- L2	L	78	59
SI-11-2063- L3	M	76	57
SI-11-2064- L4	NM	72	57
SI-11-2065- L5	S	69	50

**Stainless steel “Perforated” impression trays
with retention rim.**



Set of 10 pieces.

SI-11-2103

set
10

For Upper Jaw



		mm.	mm.
SI-11-2151- U1	XL	81	63
SI-11-2152- U2	L	74	60
SI-11-2153- U3	M	71	55
SI-11-2154- U4	NM	67	53
SI-11-2155- U5	S	65	51

For Lower Jaw



		mm.	mm.
SI-11-2161- L1	XL	80	63
SI-11-2162- L2	L	78	59
SI-11-2163- L3	M	76	57
SI-11-2164- L4	NM	72	57
SI-11-2165- L5	S	69	50

INOX
18/08

Stainless steel “Solid” impression trays with retention rim.


Impression Trays

Abdrucklöffler

Portamoldes

Porte-empreintes


Porta-impronte



set
6

Set of 6 pieces.


SI-11-2301 U4



set
8

Set of 8 pieces.

SI-11-2311




set
12

Set of 12 pieces.

SI-11-2304

For Upper Jaw



For Lower Jaw



		mm.	mm.
SI-11-2311- U1	XXS	61	47
SI-11-2312- U2	XS	64	51
SI-11-2313- U3	S	66	53
SI-11-2314- U4	M	69	57
SI-11-2315- U5	L	71	60
SI-11-2316- U6	XL	73	62

		mm.	mm.
SI-11-2331- L1	XXS	55	47
SI-11-2332- L2	XS	60	49
SI-11-2333- L3	S	72	54
SI-11-2334- L4	M	74	57
SI-11-2335- L5	L	76	59
SI-11-2336- L6	XL	78	61

Partial impression trays for crown and bridge work.



SI-11-2361- U1P



SI-11-2362- L1P




SI-11-2342- L2D
SI-11-2344- L4D
SI-11-2346- L6D



SI-11-2352- L2S
SI-11-2354- L4S
SI-11-2356- L6S




**Stainless steel “Perforated”
impression trays with retention rim.**

Impression Trays
Abdrucklöffler
Portamoldes
Porte-empreintes
Porta-impronte




Set of 6 pieces.
SI-11-2401 U1

set
6



Set of 8 pieces.
SI-11-2402


set
8



Set of 12 pieces.

set
12

For Upper Jaw



For Lower Jaw



		mm.	mm.
SI-11-2411- U1	XXS	61	47
SI-11-2412- U2	XS	64	51
SI-11-2413- U3	S	66	53
SI-11-2414- U4	M	69	57
SI-11-2415- U5	L	71	60
SI-11-2416- U6	XL	73	62

		mm.	mm.
SI-11-2431- L1	XXS	55	47
SI-11-2432- L2	XS	60	49
SI-11-2433- L3	S	72	54
SI-11-2434- L4	M	74	57
SI-11-2435- L5	L	76	59
SI-11-2436- L6	XL	78	61

**Partial impression trays for crown
and bridge work.**



SI-11-2461- U1P



SI-11-2462- L1P



SI-11-2442- L2D
SI-11-2444- L4D
SI-11-2446- L6D



SI-11-2452- L2S
SI-11-2454- L4S
SI-11-2456- L6S


Stainless Steel

INOX
18/08

Anteriorly depressed with retention rim.
Stainless Steel "Solid"

"Solid" stainless steel

For Upper Jaw



set
6

Set of 6
pieces.
SI-11-2601



		mm.	mm.
SI-11-2611- U1	L	81	67
SI-11-2614- U4	M	71	64
SI-11-2617- U7	S	70	57

For Lower Jaw





		mm.	mm.
SI-11-2631- L20	L	81	60
SI-11-2634- L21	M	79	57
SI-11-2637- L22	S	76	54

Anteriorly depressed with retention rim.
Stainless steel “Perforated”




“Perforated” stainless steel

For Upper Jaw



		mm.	mm.
SI-11-2711- U1	L	81	67
SI-11-2714- U4	M	71	64
SI-11-2717- U7	S	70	57


For Lower Jaw



		mm.	mm.
SI-11-2731- L20	L	81	60
SI-11-2734- L21	M	79	57
SI-11-2737- L22	S	76	54

INOX
18/08


Stainless steel “Solid” impression trays without retention rim.



set
6

Set of 6 pieces.


SI-11-2901 U4



set
8

Set of 8 pieces.

SI-11-2902




set
12

Set of 12 pieces.

SI-11-2904

For Upper Jaw



For Lower Jaw



		mm.	mm.
SI-11-2911- U1	XXS	61	47
SI-11-2912- U2	XS	64	51
SI-11-2913- U3	S	66	53
SI-11-2914- U4	M	69	57
SI-11-2915- U5	L	71	60
SI-11-2916- U6	XL	73	62

		mm.	mm.
SI-11-2931- L1	XXS	55	47
SI-11-2932- L2	XS	60	49
SI-11-2933- L3	S	72	54
SI-11-2934- L4	M	74	57
SI-11-2935- L5	L	76	59
SI-11-2936- L6	XL	78	61

Partial impression trays for crown and bridge work.



SI-11-2942- L2D
SI-11-2944- L4D
SI-11-2946- L6D



SI-11-2952- L2S
SI-11-2954- L4S
SI-11-2956- L6S




SURGIFUSION

INTERNATIONAL



Stainless steel “Perforated”
impression trays without
retention rim.

Impression Trays
Abdrucklöffel
Portamoldes
Porte-empreintes
Porta-impronte




Set of 6 pieces.
SI-11-3101 U1

set
6



Set of 8
pieces.
SI-11-3102


set
8



Set of 12 pieces.
SI-11-3104

set
12

For Upper Jaw



For Lower Jaw



		mm.	mm.
SI-11-3111- U1	XXS	61	47
SI-11-3112- U2	XS	64	51
SI-11-3113- U3	S	66	53
SI-11-3114- U4	M	69	57
SI-11-3115- U5	L	71	60
SI-11-3116- U6	XL	73	62

		mm.	mm.
SI-11-3131- L1	XXS	55	47
SI-11-3132- L2	XS	60	49
SI-11-3133- L3	S	72	54
SI-11-3134- L4	M	74	57
SI-11-3135- L5	L	76	59
SI-11-3136- L6	XL	78	61

Partial impression trays for crown
and bridge work.






SI-11-3142- L2D
SI-11-3144- L4D
SI-11-3146- L6D




SI-11-3152- L2S
SI-11-3154- L4S
SI-11-3156- L6S

EHRICKE Stainless steel impression trays.



For Partially Toothed Upper Jaws

	mm.	mm.
SI-11-3711- Op1	62	50
SI-11-3712- Op2	70	56



For Partially Toothed Lower Jaws

	mm.	mm.
SI-11-3731- Up1	70	55
SI-11-3732- Up2	78	59




For Edentulous Upper Jaw

	mm.	mm.
SI-11-3811- Ou1	64	46
SI-11-3814- Ou2	68	53
SI-11-3817- Ou3	79	61




For Edentulous Lower Jaw

	mm.	mm.
SI-11-3831- Uu1	70	48
SI-11-3834- Uu2	70	53
SI-11-3837- Uu3	72	57



set
8

Set of 8 pieces.
SI-11-3702



set
6


Set of 6 pieces.
SI-11-3801

CH 12


Pliers For Orthodontics
Zangen für Kieferorthopädie
Pinzas para orthodoncia
Pinces pour orthodontie
Pinze Ortodontiche



Flat Nose Pliers



Bending Pliers





SI- 12-1512
12,5 cm/ 5"




SI- 12-1511
11 cm/ 4 1/2"



SI- 12-1611
11 cm/ 4 1/2"




Diamant
SI- 12-1716
16 cm/ 6 1/4"




Stadler
SI- 12-1815
15 cm/ 6"


Wire, Arch and clasp bending pliers




Young
SI- 12-1913
13 cm/ 5"




Aderer
SI- 12-1911
11,5 cm/ 4 1/2"



Aderer
SI- 12-2013
13 cm/ 5"




Angle
SI- 12-2112
12 cm/ 4 3/4"




Adams
SI- 12-2113
13 cm/ 5"


Crown Pliers, Contouring Pliers



Reynolds
SI- 12-2213
13 cm/ 5"




Johnson
SI- 12-2313
13 cm/ 5"



Abell
SI- 12-2413



How
SI- 12-2414
14 cm/ 5 1/2"



How
SI- 12-2512 12,5 cm/ 5"
SI- 12-2514 14 cm/ 5 1/2"





SURGIFUSION

INTERNATIONAL Pliers, Copper Ring and
Bridge Removing Pliers


Pliers For Orthodontics
Zangen für Kieferorthopädie
Pinzas para orthodoncia
Pinces pour orthodontie
Pinze ortodontiche



Waldsachs-Universal
SI- 12-2615
15,5 cm/ 5"



Universal
SI- 12-2715
15,5 cm/ 5"



SI- 12-2815
15,5 cm/ 5"



Lodo
SI- 12-3015
15 cm/ 5"



Furrer
SI- 12-2915
15,5 cm/ 5"



SI- 12-3013
13 cm/ 5¼"



Lingual Arch Plier
SI- 12-3112
11.5 cm/ 4½"
Designed to form double and triple bends in either .030" (0.76 mm) or .036" (0.91 mm) wire for inserting lingual sheaths



SI- 12-3212
12 cm/ 4¾"

Remover



SI- 12-3213
13 cm/ 5"

Remover



De La Rosa
SI- 12-3312
11.5 cm/ 4½"
Arch Forming & Contouring Plier
With three grooves, for wires .016" (0.41 mm), 0.81" (0.46 mm) and .022" (0.56 mm).



De La Rosa
SI- 12-3412
11.5 cm/ 4½"
Arch Forming Plier, without grooves
Same Plier but without grooves



Remover

D.B. Bracket Removing Plier
SI- 12-3512 T/C
12 cm/ 4¾"
Useful for ceramic and stainless steel brackets. The tips wedge between edge of bracket and tooth surface and remove it with ease.



Remover

SI- 12-3612 T/C
12 cm/ 4¾"

Bond Bracket



SI- 12-3712 T/C
12 cm/ 4¾"



SI- 12-3812 T/C
12 cm/ 4¾"



SI- 12-3912 T/C
12 cm/ 4¾"




D.B. Bracket Removing Plier, with Pad
SI- 12-4012 T/C
12 cm/ 4¾"
Designed for quick and comfortable removal of direct bond brackets. It will also remove adhesive from the tooth.




Posterior Band Removing Plier, Long
SI- 12-4112 T/C
12 cm/ 4¾"
Will remove band with maximum patient comfort. The tip is postured in the middle of the pad for band removal ease. A plastic padded tip and a sharp removing beak allow, with a slight pressure, easy band removal.

Wire Bending Pliers



SI- 12-4212
12 cm/ 4¾"
Medium Wire Plier,
Round & Concave Plier
For wires up to .28" (0.71 mm)




SI- 12-4312
11.5 cm/ 4½"
Double Rounded Jaw Plier
For wires up to .030".




SI- 12-4412
12 cm/ 4¾"


Wire Bending Pliers



Jarabak
SI- 12-4413
12.5 cm/ 5"
Excellent for precise bending of light wires. Precision grooves assure accurate bending and closing loops. Wires up to .020".



Reynolds
SI- 12-4512
12 cm/ 4¾"



Light Wire Plier
SI- 12-4513
12.5 cm/ 5"
The beaks are more tapered and longer than Angle Plier. Excellent for bending small diameter loops. Used for wires up to .016" (0.41 mm)

Wire Arch Bending & Forming Pliers



Tweed
SI- 12-4611 T/C
11.5 cm/ 4½"
Rectangular Arch Pliers, Tweed Style
Forming plier for square or rectangular wires. The blades are designed to be parallel at .020" (0.51 mm) opening. The working edges are carefully diamond-honed to prevent wire scoring.




Tweed (Angle)
SI- 12-4612 T/C
12 cm/ 4¾"



Tweed (Angle)
SI- 12-4711 T/C
11 cm/ 4½"



SI- 12-4811 T/C
11 cm/ 4½"



Aderer
SI- 12-4812
11.5 cm/ 4½"
Three Jaw Thin Plier,
For wires up to .028" (0.71 mm).



Aderer
SI- 12-4912
12 cm/ 4¾"
Three Jaw Plier, Aderer
The precision aligned tips ensure for consistent bends. Gradually tapered for delicate bending. For wires up to 0.036" (0.91 mm)

Wire Arch Bending & Forming Pliers



How
SI- 12-4913
13 cm/ 5¼"
Utility Plier , How Long well tapered shape, With 3/32" (2.4mm) gripping pads which are serrated and shaped , so that their centre line coincides with the axis of rotation of the plier and holds the wire securely. Beaks are rounded to be completely lip safe.




How
SI- 12-5013
12.5 cm/ 4¾"
Utility Plier , How Style With 1/8" (3.2 mm) gripping pads, which are serrated to hold wire securely. Beaks are completely lip safe.




Weingart
SI- 12-5113 T/C
12.5 cm/ 4¾"
Utility Plier, Weingart Classic design, generally utility plier. Pointed serrated tips. Holds wire firmly at convenient angle. The beaks are rounded to make the plier completely mouth safe. Useful for removal of ceramics

Loop Forming /Bending Pliers



Adams
SI- 12-5212
12 cm/ 4¾"



Angle
SI- 12-5312
11.5 cm/ 4½"
Bird Beaks Pliers, Angle It is used for working round wires up to .030" (0.76 mm) in diameter. The beaks are .025" (0.64 mm) at the tip. The beaks are parallel at 0.20" (0.51 mm) opening.

Loop Forming /Bending Pliers



Mckellops Plier
SI- 12-5313
12.5 cm/ 5¼"
For wires up to .059" (1.5 mm).



Tweed Loop
Forming



Tweed
SI- 12-5413
12.5 cm/ 5¼"

Loop Forming, Tweed Style Pliers Designed for making loops in round or rectangular wires up to 0.022x0.28" (0.56 x 0.71 mm). The cylindrical beak has three sections: 0.045" (1.14mm), .060" (1.52mm), .075" (1.90mm). Each section is 0.090" in length. The concave beak has slight parallel serrations which orient and keep the wire at a right angle to the beaks while the loop is being formed.



Omega Loop
Forming



Tweed
SI- 12-5513
12.5 cm/ 5¼"
Omega Loop Forming Pliers, Tweed Style
It just differs from the Model 12-5413 in the concave beak which is rounded for making more uniform omega loops.



Single Diameter
Loop Forming



Tweed
SI- 12-5613
12.5 cm/ 5¼"



SI- 12-6612
Loop Closing Plier, Nance Style
Box joint with 4 sectioned tip for forming
various size loops in rectangular and round wires.
For wires up to .028" (0.71 mm).

Wire Cutting Pliers



SI- 12-5713 T/C

12.5 cm/ 5"
Hard Wire Cutter, Flat Back
Cutting capacity up to .020 round,
Arch wire up to .022 x .028 and all
sizes of ligature wire and pins.



SI- 12-5813 T/C

12.5 cm/ 5"



SI- 12-5913 T/C

13 cm/ 5½"



SI- 12-6012 T/C

12 cm/ 5"



SI- 12-6112 T/C

11.5 cm/ 4½"



SI- 12-6212 T/C

12 cm/ 5"
Mini cut & hold distal end cutter
cuts with in 0.5mm from end of
buccal tube. Use on wire up to
0.22" x 0.28" rectangular and
0.20" round.



SI- 12-6312 T/C

11.5 cm/ 4½"
Distal End Cutter with safety hold
It shear cuts hard wire and braided
wire and then, safely holds the loose
distal ends of the wire, thus preventing
from falling down into the oral tissues.
The cutting edges are diamond honed.
Maximum cutting capacity: rectangular
and twist wires up to .022 x .028"
(0.56 x 0.71); round wire up to .020"
(0.51 mm) in diameter.

CH 13-18

CH 13	Knives, Scalpels Messer, Skalpelle Bisturies Bistouris Bisturi
CH 14	Scissors Scheren Tijeras Ciseaux Forbici
CH 15	Forceps, Hemostatic Forceps Pinzetten, Arterienklemmen Pinzas, Pinzas hemostáticas Pinces, Pinces hémostatiques Pinzette, Pinze emostatiche
CH 16	Needle holders Nadel halter Portaguja Porte-aiguilles Porta aghi
CH 17	Retractors Wundhaken Separadores Ecarteurs Retrattori
CH 18	Bone Cutting Forceps, Bone Rongeurs Hohlemeißelzangen, Knochen Splitterzangen Pinzas Cortentes para huesos Pinces à os coupantes, Pinces-gouges Scalpellini e pinze alveolotome





Ash
SI-25.01
SI-25.02




1/1
SI-25.01
Fig.1



SI--25.02
Fig.2



SI-25.03
Fig.3



1/1
SI-03.10
Fig.10



SI-03.11
Fig.11



Standard
SI-04.03
NO.3



SI-09.04
NO.4



SI-06.17
NO.7K



SI-05.07
NO.7



SI-10.40
NO.3
Straight




1/1
SI-03.12
Fig.12



SI-03.13
Fig.13



1/1
SI-03.14
Fig.14



SI-03.15
Fig.15

Operating Scissors



Standard
SI-14.14
SI-14.64 "TC"
SI-18.14
SI-18.64 "TC"
14,5 cm/ 5¾"
Blunt/ Blunt



Standard
SI-15.14
SI-14.64 "TC"
SI-18.14
SI-18.64 "TC"
14,5 cm/ ¾"
Sharp/ Sharp



Standard
SI-16.14
SI-14.64 "TC"
SI-17.14
SI-17.64 "TC"
14,5 cm/ 5¾"
Sharp/ Blunt



Metzenbaum
SI-36.14
14,5 cm/ 5¾"



Metzenbaum
SI-37.14
14,5 cm/ 5¾"

Gum Scissors
Zahnfleischscheren
Tijeras para encias
Ciseaux à gencives
Forbici per gengiva



Metzenbaum-Delicate
SI-42.14 14,5 cm/ 5¾"
SI-44.18 18 cm/ 7"




Metzenbaum-Delicate
SI-43.14 14,5 cm/ 5¾"
SI-45.18 18 cm/ 7"



SI-01.11-SI-08.91
11,5 cm/ 4½"




SI-01.11



SI-02.11




SI-08.90




SI-08.91

Gum Scissors
Zahnfleischscheren
Tijeras para encias
Ciseaux à gencives
Forbici per gengiva




SI-90.13
13 cm/ 5"

SI-91.13
13 cm/ 5"



Wagner
SI-95.12
12 cm/ 4¾"

SI-96.12
12 cm/ 4¾"



Dean
SI-25.17
17 cm/ 6¾"




Quinby
SI-110.13
12,5 cm/ 5"


Quinby
SI-111.13
12,5 cm/ 5"



La Grange
SI-120.11
10,5 cm/ 4"



Locklin
SI-115.16
16 cm/ 6¼"



Universal
SI-23.12
12 cm/ 4¾"

Wire Cutting Scissors


Scissors, with tungsten carbide Inserts

Scheren mit hartmetall-Einsätzen

Tijeras con plaguitas de carburo de tungsteno

Ciseaux avec plaquettes en carbure de tungstène

forbici con placchette in carburo di tungsteno




SI-91.53
13 cm/ 5"

Locklin
SI-115.66
16 cm/ 6¼"

Kelly
SI-46.66 16 cm/ 6¼"
SI-47.66 16 cm/ 6¼"

Metzenbaum
SI- 36.64 14 cm/ 5½"
SI- 37.64 14 cm/ 5½"
SI- 38.70 20 cm/ 8"
SI- 39.70 20 cm/ 8"
Metzenbaum-Fine
SI- 42.64 14 cm/ 5½"
SI- 43.64 14 cm/ 5½"
SI- 44.70 20 cm/ 8"
SI- 45.70 20 cm/ 8"



SI- 08.92
11,5 cm/ 4½"

Ragnell (Killner)
SI- 59.63
13 cm/ 5½"

Universal
SI- 24.12
12 cm/ 4¾"

Dean
SI- 25.57
17 cm/ 6¾"

Locklin
SI- 115.66
16 cm/ 6¼"
Saw Edge



	Standard	Standard	Adson
12 cm/ 4¾"	DS-01.11	DS-05.11	DS-15.12
14.5cm/ 5¾"	DS-01.14	DS-05.14	
15cm/ 6"	DS-01.15	DS-05.15	DS-15.15
16 cm/ 6¼"	DS-01.16	DS-05.16	
18 cm/ 7"	DS-01.18	DS-05.18	
20 cm/ 8"	DS-01.20	DS-05.20	
25 cm/ 10"	DS-01.25	DS-05.25	



Semken
DS-30.15
15.5 cm/ 6"



Adson
DS-16.12 12 cm/ 4¾"
DS-16.15 15 cm/ 6"




Adson Brown
DS-21.12 12 cm/ 4¾"
DS-21.15 15 cm/ 6"

With
Tungsten Carbide
"TC"



	Standard "TC"	Narrow "TC"	Adson "TC"
12 cm/ 4¾"	DS-01.62	DS-03.62	DS-15.62
14.5cm/ 5¾"	DS-01.64	DS-03.64	
15cm/ 6"	DS-01.65		DS-15.65
16 cm/ 6¼"	DS-01.66	DS-03.66	
18 cm/ 7"	DS-01.68	DS-03.68	

Cushing "TC"
DS-52.67
17cm / 6¾"



London College
15 cm/ 6"
SI- 15-0915 -SI-15-1515
WO: Without Lock
SI- 15-0916 -DS-15-1516
WL: With Lock

Flagg
16 cm/ 6¼"
SI- 15-1616 Smooth
SI- 15-1626 Serrated



Meriam
16 cm/ 6¼"
SI-15-1716 - SI-15-2216
WO: Without Lock
SI- 15-1717 -SI-15-2217
WL: With Lock

Perry
13 cm/ 5"
SI- 15-2410 Serrated
SI- 15-2413 Smooth



Albrecht
SI- 15-2516 Sharp
SI- 15-2517 Blunt
16 cm/ 6 1/4"



USA-Mod.
SI- 15-2615
15 cm/ 6"



SI- 15-2715
15 cm/ 6"



Caulk
SI- 15-2815
16,5cm/ 6 1/2"



SI- 15-2915
15 cm/ 6"

Articulating Forceps



Miller
SI- 15-3415
15cm/ 6 "



USA-Mod.
SI- 15-3016
16 cm/ 6 1/4"



USA-Mod.
SI- 15-3116
16 cm/ 6 1/4"



USA-Mod.
SI- 15-3216
16 cm/ 6 1/4"



Crane-Keplan
SI- 15-3301
14 cm/ 5 1/2"



Crane-Keplan
SI- 15-3302
14 cm/ 5 1/2"


Hemostatic Forceps

Arterienklemmen

Pinzas hemostaticas

Pinces hémostatiques

Pinze emostatiche



1/1
SI- 40.09


SI- 41.09

Hartmann

SI- 40.09 - SI- 41.09

10 cm/4"
Straight

10 cm/4"
Curved



SI- 10.12


SI- 11.12

Halstead

SI- 10.12 - SI- 11.12

12 cm/4 3/4"
Straight

12 cm/4 3/4"
Curved



1 : 2
SI- 18.12

1 : 2
SI- 19.12

Halsted-Mosquito

SI- 18.12 - SI- 19.12

12,5 cm/5"
Straight

12 cm/5"
Curved



Jones

SI- 19.06 6 cm/ 2 1/2"

SI- 20.09 9 cm/ 3 1/2"




Backhaus

SI- 25.08 8 cm/ 3 1/4"

SI- 25.10 11 cm/ 4 1/2"

SI- 26.13 13 cm/ 5"



SI- 11.20
SI- 13.26

SI- 12.20
SI- 14.26


Gross-Maier

SI- 11.20 - SI- 13.26

SI- 12.20 - SI- 14.26

20cm/ 8"

26cm/ 10 1/4"





Frazier
SI- 11.12-
SI- 11.13
13 cm/ 5 1/4"



Kilner
SI- 11.15
15 cm/ 6"




Guthrie
SI- 25.13
13 cm/ 5"



Gillies
SI- 12.19-
SI- 12.20
19 cm/ 7"



Senn-Miller
SI- 01.16
16 cm/ 6 1/4"




SI- 31.01 -
SI- 31.04
16 cm/ 6 1/4"



Farabeuf
SI- 46.12
12 cm/ 4 3/4"



Farabeuf
SI- 47.15
15 cm/ 5"



Collin-Hartmann
SI- 45.14 14 cm/ 5 1/2"
SI- 45.16 16 cm/ 6 1/4"



Roux
Roux Set of 3.
SI- 49.00 Set of 3
SI- 49.01 Fig.1
SI- 49.02 Fig.2
SI- 49.03 Fig.3





Fig.1 = 24 x 20 mm
Fig.2 = 29 x 26 mm
Fig.3 = 33 x 29 mm



Fig.1 = 30 x 28 mm
Fig.2 = 36 x 36 mm
Fig.3 = 38 x 43 mm



Volkmann
DS-57.21-RS-57.24
21.5 cm/ 8 ½"



Langenbeck
SI- 71.01-DS-71.03
21cm/ 8 ¼ "



Obwegeser
DS-12.40-
DS-12.80
22 cm/8 ¾ "

Obwegeser
DS-10.05 -
DS-10.07



Williger
SI- 11.17
17 cm/6 ¾"



Williger
SI- 12.04
16 cm/6 ¼"




Freer
SI- 28.21
21.5 cm/ 8 ½"



SI- 19.20
17,5 cm/7"



Molt
SI- 19.21
19 cm/7 ½"




Lengenbeck
SI- 37.35
19 cm/7 ½"

Bone Nibblers.



1mm
SI- 37.01
2mm
SI- 37.02
3mm
SI- 37.03

Citelli
**SI- 37.01-
SI- 37.03**
17 cm/6 ¾"



4 x 4 mm
SI- 36.01

5 x 5 mm
SI- 36.02



6 x 6 mm
SI- 36.03



Kerrison
**SI- 36.01-
SI- 36.03**



Friedman
SI- 01.14
14 cm/5 ½"



- 2 mm ○ SI- 10.62
- 3 mm ○ SI- 10.63
- 4 mm ○ SI- 10.64
- 5 mm ○ SI- 10.65
- 6 mm ○ SI- 10.66
- 7 mm ○ SI- 10.67
- 8 mm ○ SI- 10.68



Keyes
SI- 10.62-
SI- 10.68




Farrell Ø 0.9 mm
SI- 16.12 12 cm/ 4 ¾"
SI- 16.14 14 cm/ 5 ½"
SI- 16.16 16 cm/ 6 ¼"




Bowmann
DS-94.00
Set of 4
SI- 94.01 Fig 1+2
SI- 94.03 Fig 3+4
SI- 94.05 Fig 5+6
SI- 94.07 Fig 7+8



Ø 1.5 mm	Ø 2 mm	
SI- 01.10	SI- 02.10	11 cm/ 4 ¼"
SI- 01.13	SI- 02.13	13 cm/ 5"
SI- 01.14	SI- 02.14	14 cm/ 5 ½"
SI- 01.16	SI- 02.16	16 cm/ 6 ¼"




Luniatchek
SI- 30.19
19 cm/7 ½"



Luniatchek
SI- 30.20
19 cm/7 ½"



Sternberg
SI- 03.15
14 cm/5 ½"



Sternberg
SI- 03.16
16 cm/6 ¼"



Bruenings
SI- 05.18
19 cm/7 ½"




SI- 03.16
Wood
16 cm/ 6 ¼"




Kroenig
SI- 03.14
15 cm/6"



Mayo
SI- 04.17
17 cm/6 ¼"




Oringer
SI- 01.00 Fig.1
SI- 01.02 Fig.2
SI- 01.03 Fig.3
SI- 01.04 Fig.4
SI- 01.13 Set of 4



Vestibulum-Haken
SI- 01.15
15 cm / 6"



Heister
SI- 05.05
13 cm / 5"



Molt
SI- 07.15
15 cm / 6"




Jennings
SI- 09.09 9 cm / 3 1/2"
SI- 09.11 11 cm / 4 1/4"
SI- 09.13 13 cm / 5"
SI- 09.15 15 cm / 6"



Roser-König
SI- 13.45 16 cm / 6 1/4"
SI- 13.50 20 cm / 7 1/2"



Denhart
SI- 14.60
13 cm / 5"



Passow
SI- 37.11-
SI- 37.14
 12 cm/ 4 3/4 "




Partsch
 13,5 cm/5 1/2 " 17 cm/6 3/4 "

	13,5 cm/5 1/2 "	17 cm/6 3/4 "	13,5 cm/5 1/2 "	17 cm/6 3/4 "
3 mm	SI- 05.03	SI- 05.53	SI- 06.03	SI- 06.53
4 mm	SI- 05.04	SI- 05.54	SI- 06.04	SI- 06.54
5 mm	SI- 05.05	SI- 05.55	SI- 06.05	SI- 06.55
6 mm	SI- 05.06	SI- 05.56	SI- 06.06	SI- 06.56



Partsch
 13,5 cm/5 1/2 " 17 cm/6 3/4 "




Cottle
SI- 28.04
SI- 28.12

SI- 14.24 Fig.0000
SI- 14.33 Fig.000
SI- 14.42 Fig.00
DS-14.50 Fig.0
SI- 14.51 Fig.1
SI- 14.52 Fig.2
SI- 14.53 Fig.3
SI- 14.54 Fig.4
SI- 14.55 Fig.5
SI- 14.56 Fig.6



Volkman
SI- 14.22- SI- 14.56
 17 cm/6 3/4 "




Williger

SI- 21.10
SI- 21.11
SI- 21.12



Partsch
SI- 30.15
 15 cm/6 " Fig .1



Partsch
SI- 31.15
 15 cm/6 " Fig .2


Bone Curettes and Mallets

Scharte Iffel und Hämmer

Curetas(Cucharillas) Cortantes y Martillos y Mazos

Curettes tranchantes et maillets

Cucchiaini per osso e martelli




SI-18-6200 - SI-18-6211
SI-18-6311 - SI-18-6333



SI- 01.16
16,5 cm/ 6 1/2"



Collin
SI- 05.20
20 cm/ 8"



Mead
SI- 06.19
Massive - solid
17 cm/ 6 3/4"

Bone Currettes

Scharte Iffel

Curetas(Cucharillas) Cortantes


Curettes tranchantes

Cucchiaini per osso



1/2

Miller-Colburn
SI-18-6411 -
SI-18-6444



1/1


SI- 18-6411

SI- 18-6422

SI- 18-6433

SI-18-6444

Miller-Colburn



1/1

SI- 18-6501

SI- 18-6502

SI- 18-6503

SI- 18-6504

Lucas

SI- 18-6401

Lucas
SI-18-6501 -SI-18-6504


Tongue Depressors / Mouth Gags / Tongue Holding forceps

Zungenspatel / Mundsperrer / Zungenzengen


Bajaleguas / Abrebocas / Pinzas ira lengne

Ecarteurs des levres / Ouvre-bouches / Pinces tire langue

Spatola per lingua / Apribocca / Pinza tira lingua



Wieder
SI- 02.21-
SI- 02.22



Needle Cases



SI- 15.05 = 50 x 25 x 4,5 mm
SI- 15.08 = 65 x 42 x 8 mm
SI- 15.06 = 75 x 25 x 4,5 mm





SI- 16.05 = 50 x 25 x 4,5 mm
SI- 16.08 = 65 x 42 x 8 mm
SI- 16.06 = 75 x 25 x 4,5 mm


Trays





SI- 05.63
235 x 190 x 40 mm



complete	insert only
SI- 19.65 Ø = 65 mm h = 15 mm	SI- 19.70 Ø = 62 mm h = 13 mm



SI- 05.82
Ø 200 x 70 mm



	l	b	h
SI- 04.20	200	150	38 mm
SI- 04.25	250	200	50 mm
SI- 04.30	300	250	50 mm
SI- 04.35	355	305	50 mm
SI- 04.40	405	305	50 mm

SI- 04.20- SI- 04.40

Forceps Jars



Size in mm Ø X h	Jar without lid	Jar with lid
30 x 90	SI- 15.30	SI- 15.54
50 x 100	SI- 15.50	SI- 15.55
50 x 130	SI- 15.51	SI- 15.56
50 x 175	SI- 15.52	SI- 15.57
75 x 100	SI- 15.75	SI- 15.80
75 x 130	SI- 15.76	SI- 15.81
75 x 175	SI- 15.77	SI- 15.82
100 x 175	SI- 15.98	SI- 15.88
100 x 200	SI- 15.99	SI- 15.89

Instrument Box



SI- 20.20
205 x 105 x 40 mm

Kidney Trays



Size in mm l x b x h		
170 x 100 x 35	SI- 10.17	17 cm/ 6 ¾"
250 x 140 x 40	SI- 10.25	25 cm/ 10"
275 x 150 x 45	SI- 10.27	27 cm/ 10 ¾"

Instrument Boxes



Size in mm
L x w x h

165 x 85 x 35	SI- 21.46
200 x 100 x 35	SI- 21.50
200 x 100 x 50	SI- 21.51
230 x 130 x 50	SI- 21.53
260 x 150 x 50	SI- 21.56

270 x 235 x 50	SI- 21.57
300 x 200 x 50	SI- 21.60
325 x 275 x 30	SI- 21.62
325 x 275 x 50	SI- 21.63
420 x 80 x 40	SI- 21.72

420 x 175 x 50	SI- 21.75
----------------	------------------

Cotton Dispensers



Methot
SI- 35.01 **SI- 35.02**
ø50 x 75 mm
With lid Without lid



Spirit Lamp
SI- 19-1505
120 cc

Cotton wool containers



Size in mm Lid with knob

Ø X h	
75 x 60	SI- 08.07
75 x 75	SI- 08.08
100 x 60	SI- 08.10
100 x 100	SI- 08.11
125 x 80	SI- 08.12
125 x 125	SI- 08.13
150 x 100	SI- 08.14
150 x 150	SI- 08.15
200 x 100	SI- 08.20
200 x 150	SI- 08.21
200 x 200	SI- 08.22
240 x 160	SI- 08.24
240 x 230	SI- 08.25
290 x 160	SI- 08.29
290 x 240	SI- 08.30

Sterilization Drums

Sterilization, Asepsis

Sterilisation, Asepsis

Esterilización, Asepsia

Stérilisation, Asepsis

Sterilizzazione, Asepsi




Schimmelbusch
(Horizontal Ventilation)

Ø	h	
140 x 140 mm		SI- 25.14
160 x 150 mm		SI- 25.16
190 x 160 mm		SI- 25.19
190 x 190 mm		SI- 25.20
240 x 120 mm		SI- 25.24
240 x 160 mm		SI- 25.25
240 x 240 mm		SI- 25.27
290 x 160 mm		SI- 25.29
290 x 240 mm		SI- 25.30
340 x 160 mm		SI- 25.34
340 x 240 mm		SI- 25.35
390 x 140 mm		SI- 25.39
390 x 190 mm		SI- 25.40
390 x 290 mm		SI- 25.41



Konrich
(Vertical Ventilation)

Ø	h	
140 x 140 mm		SI- 26.14
160 x 150 mm		SI- 26.16
190 x 160 mm		SI- 26.19
190 x 190 mm		SI- 26.20
240 x 120 mm		SI- 26.24
240 x 160 mm		SI- 26.25
240 x 240 mm		SI- 26.27
290 x 160 mm		SI- 26.29
290 x 240 mm		SI- 26.30
340 x 160 mm		SI- 26.34
340 x 240 mm		SI- 26.35
390 x 140 mm		SI- 26.39
390 x 190 mm		SI- 26.40
390 x 290 mm		SI- 26.41




Norm Trays


SI-19-2101 - SI-19-2104
288 x 187 x 40 mm

SI-19-2201 - SI-19-2204
288 x 93 x 40 mm

Inserts for Trays



SI-19-2321



SI-19-2322

SI-19-2323



SI-19-2422




SI-19-2423

CH 20

Surgery
Chirurgie
Chirgula
Chirurgie
Chirurgia



With additional side holes



6 mm
round

6 mm
oval



Nubor
SI- 20-0121
21 cm/8½"



Nubor
SI- 20-0122
21 cm/8½"



Coupland
SI- 20-0217
17 cm/ 6 ¾"



- SI- 20-0301
1,8 mm
- SI- 20-0302
2,8 mm
- SI- 20-0303
3,8 mm

Byrd
SI- 20-0301-
SI- 20-0303
20 cm/ 8"



Cogswell
SI- 20-0423
23 cm/ 9"



Adson
SI- 39.32
4 mm
21 cm/ 8 ½"



Frazier
Charr. / Fr.
SI- 36.86 6
SI- 36.88 8
SI- 36.90 10
SI- 36.92 12
19 cm/ 7 ½"



Yankauer
SI- 30.00
27 cm/ 11"



SI- 20-0702-1.8cc
SI- 20-0703-2.2cc
(Complete with 2 tips)



SI- 20-0801



SI- 20-0802



SI- 20-0803



SI- 20-0804




SI- 20-0805



SI- 20-0806

Chiropody Instruments



Sterilization Forceps



Spatulas



Cheatle
SI- 03.27
27 cm / 10 1/2"

SI- 19.30
30 cm / 12"


SI- 20-1919 = 19 cm / 7 1/2"
SI- 20-1922 = 22 cm / 8 1/2"
SI- 20-1925 = 25cm / 10"

SI- 20-2021
21 cm / 8 1/4"
105 x 30 mm
elastic

SI- 20-2120
20 cm / 8"
100 x 23 mm
rigid

SI- 20-2219
19 cm / 7 1/2"
100 x 17 mm
elastic

Napkin Holders



Professional
SI- 20-2301
length adjustable



Adams
SI- 20-2302

Finger Guards



SI- 20-2401



SI- 20-2402



SI- 20-2403

Langenbeck
SI-20-2401 - SI-20-2403

Sand glass



SI- 20-2500

Spring Caliper




SI- 20-2601



SI- 20-2602

Micro meters



Boley
SI- 20-2710



Iwanson
SI- 20-2603



Iwanson
SI- 20-2604

Oral wedge



Maunder
SI- 64.10

0 - 100 mm

Mouth mirrors and handles for mouth mirrors

Mundspiegel und Mundspiegelgriffe

Espejos de boca y mango para espejos de boca

Miroirs à bouche Manche pour miroirs

Specchietti e manici per specchietti



Ø mm	Sizes	Plane, Cone-Socket	Magnifying, Cone Socket	Ø mm	Sizes	Plane, simple thread	Magnifying, simple thread
22	4	SI- 20-2922	SI- 20-3022	22	4	SI- 20-3122	SI- 20-3222
24	5	SI- 20-2924	SI- 20-3024	24	5	SI- 20-3124	SI- 20-3224

Handles for mirrors



1/1



Secco



With simple thread **SI- 20-3301**

With Cone-Socket **SI- 20-3302**

SI- 20-3401

SI- 20-3501

SI- 20-3601


SI- 20-3701

CH 21

Dental instruments sets
Satz Instrumente
Sistema de Instrumentos
Ensemble d' instruments
Insieme degli strumenti



Supplied in Leather wallets





SI- 21-0108, Set of 8., Containing

- SI- 01-0101** Dental Extracting Forceps, No.1
- SI- 01-0113** Dental Extracting Forceps, No.13
- SI- 01-0117** Dental Extracting Forceps, No.17
- SI- 01-0118** Dental Extracting Forceps, No.18
- SI- 01-0751** Dental Extracting Forceps, No.51A
- SI- 01-1067** Dental Extracting Forceps, No.67
- SI- 01-1379** Dental Extracting Forceps, No.39

SI- 21-0209 Set of 9, Containing

- SI- 01-0333** Dental Extracting Forceps, No.33
- SI- 01-0751** Dental Extracting Forceps, No.51A
- SI- 01-1067** Dental Extracting Forceps, No.67
- SI- 01-1379** Dental Extracting Forceps, No.79
- SI- 05-0203** Root Elevators, 3mm, 13.5 cm/ 5¼"
- SI- 05-2302** Root Elevators, No. 302
- SI- 05-2303** Root Elevators, No. 303
- SI- PJ.69.14** Mayo Hegar Needle Holder, 16 cm, T/C
- SI- PJ.23.56** Mathieu Needle Holder, 14cm, T/C



SI- 21-0309 Set of 9 , Containing


- SI- 01-0333** Dental Extracting Forceps, No.33
- SI- 01-0751** Dental Extracting Forceps, No.51
- SI- 01-1067** Dental Extracting Forceps, No.67
- SI- 01-1379** Dental Extracting Forceps, No.79
- SI- 05-0203** Root Elevators, 3mm, 13.5 cm/ 5¼"
- SI- 05-2302** Root Elevators, No. 302
- SI- 05-2303** Root Elevators, No. 303
- SI- 05-4320** Dental Extracting Forceps, No.320
- SI- 05-4321** Dental Extracting Forceps, No.321

SI- 21-0408 Set of 8, Containing

- SI- 01-0101** Dental Extracting Forceps, No.1
- SI- 01-0104** Dental Extracting Forceps, No.4
- SI- 01-0107** Dental Extracting Forceps, No.7
- SI- 01-0113** Dental Extracting Forceps, No.13
- SI- 01-0117** Dental Extracting Forceps, No.17
- SI- 01-0118** Dental Extracting Forceps, No.18
- SI- 01-0220** Dental Extracting Forceps, No.20
- SI- 01-1686** Dental Extracting Forceps, No.6C



Periodontal Instruments Set



SI- 21-0508 Set of 16 Periodontia Instruments, Containing in artificial Leather case :

- SI- 08-2101** Goldman-Fox Curettes 16cm/ 6 1/4"
- SI- 08-2102** Goldman-Fox Curettes 16cm/ 6 1/4"
- SI- 08-2103** Goldman-Fox Curettes 16cm/ 6 1/4"
- SI- 08-2104** Goldman-Fox Curettes 16cm/ 6 1/4"
- SI- 08-2105** Goldman-Fox Curettes 16cm/ 6 1/4"
- SI- 08-2106** Goldman-Fox Curettes 16cm/ 6 1/4"
- SI- 08-1507** Goldman-Fox Knives
- SI- 08-1508** Goldman-Fox Knives
- SI- 08-1509** Goldman-Fox Knives
- SI- 08-1510** Goldman-Fox Knives
- SI- 08-1511** Goldman-Fox Knives
- SI- 08-1521** Goldman-Fox Knives
- SI- PC.90.13** Goldman Fox, Gum Scissors 13 cm/ 5 1/8"
- SI- 15-3301** Crane-Keplan Dental Tweezers 14 cm/ 5 1/2"
- SI- 15-3302** Crane-Keplan Dental Tweezers 14 cm/.5 1/2"

	<i>Page</i>		<i>Page</i>
A		Crane-Kaplan Dental Tweezers	104
Abell Pliers, Contouring	89	Crane-Kaplan Gingivectomy Knives	51
Adams Napkin Holders	129	Crile-Wood Needle Holders	107
Adams Pliers, Wire, Arch and Clasp Bending	89,94	Crile-Wood Needle Holders, TC	108
Aderer Pliers, Wire, Arch and Clasp Bending	89,93	Cross Periodontia Instruments	53
Adson Forceps, Dissecting TC	102	Crown Instruments	71,72
Adson Forceps, Dressing	102	Cryer Extracting Forceps, American Pattern	23
Adson Forceps, Tissue	102	Cryer Root Elevators	28
Adson Suction Tubes	126	Cups, Finger	36
Adson-Brown Forceps	102	Curettes, Bone	115-117
Ainsworth Rubber Dam Instruments	74	Curettes, Universal	47,48
Albrecht Dental Tweezers	104	Cushing Forceps, Dissecting, TC	102
Allis Forceps, Tissue and Organ Holding	106	Cutting Instruments	66
Allis-Baby Forceps, Tissue and Organ Holding	106		
Amalgam Instruments	35-37	D	
Andrew Tongue Depressors	118	Darby-Perry Excavators	60
Angle Pliers, Wire, Arch and Clasp Bending	89,94	Darby-Perry Scalers	46
Apical Flohr Root Elevators	28	Davis-Boyle Mouth Gags	118
Apical Root Elevators	28	De La Rosa, Pliers	91
Articulating Forceps	104	Dean Gum Scissors	100
Ash Gum Knives	98	Dean Gum Scissors, TC	101
Ash Plastic Filling Instruments	39	Decimal Spring Calipers	129
Axial Filling Instruments	40,41	Denhart Mouth Gag	114
		Derf Needle Holders, TC	108
B		Diamant End and Side Nipper	89
Babcock Forceps, Tissue and Organ Holding	106	Dispenser, Gauze and Cotton	122
Baby-Crile Needle Holders	108	Dividing Compasses	129
Backhaus Towel Forceps	105	Dressing Forceps	102,105
Barnhart Universal Curettes	48	Drums, Sterilization	123
Bates Universal Curettes	48	Dumont Pliers	88
Beale Wax, Porcelain and Modelling Carvers	70		
Bein Root Elevators	27	E	
Bennet Plastic Filling Instruments	38	Ehrlicke Stainless Steel Impression Trays	85-86
Beyer Forceps, Bone Rongeur	112	Elevators & Raspatories	111
Black Cutting Instruments	66	Elevators, Root	27-32
Blumenthal Forceps, Bone Rongeur	112	Ermer Wax Knife	70
Bohler Bone Rongeurs	112	Excavators	58-61
Bohn Crown Instrument	72	Explorers	62-63
Boley Sliding Caliper	129	Extracting Forceps for Children, American Pattern	15
Bone Curettes	115-117	Extracting Forceps for Children, English Pattern	12-14
Bow Compasses	129	Extracting Forceps	2-24
Bowls, Kidney	121	Extracting Forceps, with Alfa handle	11
Bowls, Round	120	Extracting Forceps, American Pattern	16-24
Bowmann Lacrimal Probes	113	Extracting Forceps, English Pattern	2-10
Boxes, Instrument	121,122		
Brewer Rubber Dam Instruments	74	F	
Bruenings Tongue Depressors	113	F.Barry Root Elevators	30
Bunting Universal Curettes	48	Farabeuf Retractors, Soft Tissue	109
Burnishers	42	Farrell Cotton Applicators	113
Byrd Suction Tubes	126	Fedi Periodontia Instruments	55
		Fehnenstock, Wax Knife	70
C		Files	66
Calipers	129	Finger Protectors	129
Cartridge Syringes	127	Flagg Dental Tweezers	103
Cases, Needle Sterilizing	120	Flohr Root Elevators	28
Castroviejo Needle Holder	107	Forceps	102-106
Caulk Dental Tweezers	104	Forceps, Artery	105,106
Cement Spatula	43,44	Forceps, Bone Cutting	112
Cheatle Forceps, Sterilizing	128	Forceps, Bone Rongeur	112
Chip Blowers	127	Forceps, Dissecting, TC	102
Chiroprody Instruments	128	Forceps, Splinter	106
Chisels, Bone	115	Forceps, Sterilizing	128
Chompret Root Elevators	29	Forceps, Tissue and Organ Holding	106
Citelli Bone Rongeurs	111	Forceps, Tissue	102,106
Cleveland Forceps, Bone Cutting	112	Forceps, Tongue Holding	118
Cogswell Suction Tubes	126	Forceps, Towel	105
Coleman Root Elevators	29	Frahm Wax, Porcelain & Modelling Carver	68
Collin Forceps, Tongue Holding	118	Frazier Skin Hook	109
Collin Mallets	116	Frazier Suction Tubes	126
Collin-Hartmann Retractors, Soft Tissue	109	Freer Elevator	111
Columbia Scalers & Curettes	47	Friedman Forceps, Bone Rongeur	111
Containers for Cotton Wool	122	Friedman Root Elevators	26
Cottle Bone chisels	115	Furrer pliers, Copper Ring and Bridge Removing	90
Cotton Applicators	113		
Cotton Dispenser	122		
Coupland Root Elevators	27		
Coupland Suction Tubes	126		
Crane Pick Root Elevators	28		

	<i>Page</i>		<i>Page</i>
G		Lecluse Root Elevators	31
Gauze Packers	113	Locklin Gum Scissors	100,101
Gillies Needle Holder	107	Locklin Gum Scissors, TC	101
Gillies Needle Holders, TC	108	Lodo Pliers, Copper Ring and Bridge Removing	90
Gillies Skin Hooks	109	London College Tweezers, Cotton	103
Gingivectomy Knives	51	Lowrence Read, Extracting Forceps	08
Gisy-Wild Wax, Porcelain and Modelling Carver	70	Lucas Bone Curettes	117
Goldman Fox Periodontia Instruments	53,54	Luer Forceps, Bone Rongeur	112
Goldman Universal Curettes	48	Luks Nerve Instruments	65
Gracey Pereiodontia Instruments	56	Luniathek Gauze Packers	113
Gross-Maier Dressing & Cotton Swab Forceps	105		
Gum Knives	98	M	
Gum Scissors	100	Mallery Universal Curettes	48
Guthrie Skin Hooks	109	Mallets	116
Guy Tongue Holding Forceps	118	Mathieu Needle Holders	107
		Mathieu Needle Holders, TC	108
H		Matrix Bands	73
Halsey Needle Holders, TC	108	Matrix Retainers	73,74
Halsted-Mosquito Artery Forceps	105	Maunder Wedges, Oral	129
Handle for Mouth Mirrors	130	Mayo Hegar Needle Holders	107
Harris Extracting Forceps, American Pattern	17-18	Mayo Hegar Needle Holders, TC	108
Hartmann Forceps, Artery	105	Mayo Tongue Depressors	113
Heidebrink Root Elevators	29	McCall Periodontia Instruments	52
Heidemann Plastic Filling Instruments	39	McCall Scalars	46
Heister Mouth Gags	114	Mclvor Mouth Gags	118
Hemingway Mod Bone Curettes	116	Mckellops Pliers	95
Henahan Extracting Forceps, American Pattern	16	Mead Extracting Forceps	10
Hirschfeld Files	66	Mead Mallets	116
Hollenback Amalgam Instruments	37	Measuring Instruments	129
Hollenback Wax, Porcelain & Modelling Carvers	68	Meriam Tweezer, Cotton and Dressing	103
Hooks, Skin	109	Mershon Instruments for Crown Work	71
How Pliers	89,94	Methot Dispensers, Gauze and Cotton	122
Hull Extracting Forceps, American Pattern	21	Metzenbaum Delicate Scissors, Dissecting	99
Hygenist Periodontia Instruments	52	Metzenbaum Delicate Scissors, TC	101
Hylin Root Elevators	27	Metzenbaum Scissors, Dissecting	99
Hylin Wax Knife	69	Micrometers	129
		Micro-Ryder Needle Holder TC	108
I		Miller Paper Clip	71
Impression Trays	76-86	Miller Root Elevators	28
Indiana University Universal Curettes	48	Miller-Colbrun Bone Curettes	117
Instruments for Crown Work	71,72	Mirrors, Mouth	130
Interproximal Wax, Porcelain & Modelling Carvers	68	Mod. PHF Bone Curettes	116
Ivory Amalgam Instruments	36	Mod. Plastic Filling Instruments	38,40
Ivory Matrix Retainers	73,74	Mod. USA Amalgam Instruments	37
Iwanson Decimal Spring Calipers	129	Molt Mouth Gag	114
		Molt Raspatories, Periosteal	111
J		Mouth Gags	114,118
Jaquette Periodontia Instruments	52		
Jaquette Scalars	46	N	
Jarabak Pliers	93	Nabers Pernodontia Instruments	53
Jars, Forceps	121	Napkin Holders	129
Jennings Mouth Gags	114	Needle Cases	120
Johnson Pliers, Contouring	89	Needle Holders	107,108
Jones Towel Forceps	105	Needle Holders, TC	108
		Nerve Instruments	64,65
K		Nevius Extracting Forceps, American Pattern	21,22
Kells Extracting Forceps, American Pattern	22	Nippers, Cuticle	128
Kelly Gum Scissors, TC	101	Nippers, Nail	128
Kerrison Bone Rongeurs	111	Norm Trays	124
Keyes Dermal Punch	113	Nuboer Suction Tubes	126
Kidney Trays	121		
Kilner Skin Hook	109	O	
Kirkland Gingivectomy Knives	51	Obwegeser Retractors, Soft Tissue	110
Klein Extracting Forceps for Children, with Spring	14	Ochsenbein Periodontia Instruments	55
Knives & Scalpels	98	Ochsner-Kocher, Hemostatic Forceps	106
Konrich Sterilizing Drums	123	Olsen Hegar Needle Holders	107
Krallengnber Elevators	31	Olsen-Hegar Needle Holders, TC	108
Kroenig Tongue Depressors	113	Oringer Self retaining Cheek Lip Retractors	114
L		P	
La Grange Gum Scissors	100	Parmly Extracting Forceps, American Pattern	19
Langenbeck Elevators	111	Partsch Chisels, Bone	115
Langenbeck Finger Protectors	129	Passow Bone Chisels	115
Langenbeck Retractors	110	Periodontia Instruments	50-56
Langer Scalars & Curettes	47	Perry Dental Tweezers	103
Le Cron Wax, Procelain & Modelling Carvers	70	Physick Extracting Forceps, American Pattern	16

	<i>Page</i>		<i>Page</i>
Plastic Filling Instruments	38-41	Suction Tubes	126
Pliers, Bond Bracket	91,92	Surgery Instruments	125-130
Pliers, Contouring	89	Syringes Cartridges	127
Pliers, Copper Ring and Bridge Removing	90	Syringes, Water	127
Pliers, Crown	89		
Pliers, Flat Nose	88	T	
Pliers, Loop Forming & Bending	94,95	Thomas Modelling Carvers	69
Pliers, Orthodontic	88-96	Tobold Tongue Depressors	118
Pliers, Bar Bending	88	Tofflemire Matrix Bands	73
Pliers, Band Seating	91	Tofflemire-Universal Matrix Retainers	73
Pliers, Universal	90	Tomes Extracting Forceps	20
Pliers, Wire Arch and Clasp Bending	89,91-93	Tongue Depressors	113,118
Pliers, Wire Cutting Pliers	96	Towner Jaquette Periodontia Instruments	52
Pott Root Elevators	31	Trays, Half	124
Prichard Periodontia Instruments	53	Trays, Instruments Storing	120
Probes (Explorers)	62-63	Trays, Norm	124
Probes Periodontia	50	Trays, Universal	100
Professional Napkin Holders	129	Treyemann Instruments for Crown Work	71
		Tuft Wax, Porcelain and Modelling Carvers	68
Q		Turnbull Nipper, Nail	128
Quinby Gum Scissors	100	Tweed Bending & Forming Pliers	93,95
		Tweezers, Cotton & Dressing	103,104
R		U	
Ragnell (Kilner) Gum Scissors, TC	101	Universal Curettes	47,48
Ralk Forceps, Splinter	106	Universal Pliers	90
Raspatories, Periosteal	111	Universal Wire and Crown Scissors	100,101
Retractors	109-110	USA-Mod Dental Tweezers	104
Retractors, Cheek and Lip	113		
Retractors, Soft Tissue	110	V	
Reynolds Pliers	89,93	Vestibulum-Haken	114
Rochester-Pean Hemostatic Forceps	106	Volkman Bone Curettes	115
Root Elevators	26-32	Volkman Retractors	110
Roser Konig Mouth Gag	114		
Roux Retractors	109	W	
Rubber Dam Instruments	74	Wagner Gum Scissors	100
Rule Universal Curettes	48	Waldsachs-Universal Pliers	90
		Walter F.Barry Root Elevators	28,30,32
S		Ward Wax, Porcelain and Modelling Carver	68
Sachs Instruments for Crown Work	71	Ward-Kirkland Periodontia Instruments	51
Sand Glass	129	Warwick James Root Elevators	32
Scalers	46-48	Wax, Porcelain and Modelling Carvers	68-70
Scalpel Handles	98	Wedges, Oral	129
Scalpel Blades	98	Weingart Pliers	94
Schimmelbusch Sterilizing Drums	123	West Root Elevators	32
Schroeder Instruments for Crown Work	71	White Burnishers	42
Schwarz Pliers, Bar Bending	88	White Excavators	59
Scissors	99-101	White Plastic Filling Instruments	39
Scissors, Cuticle	128	Wieder Tongue Depressors	118
Scissors, Dissecting	99,101	Wiland Wax, Porcelain and Modelling Carvers	68
Scissors, Gum	99,100	Williams Plastic Filling Instruments	38
Scissors, Gum, TC	101	Williger Bone Curettes	115
Scissors, Nail	128	Williger Elevators	111
Scissors, Operating Fine	99	Winter Root Elevators	30
Scissors, Operating	99	Woodson Plastic Filling Instruments	39
Scissors, TC	99,101	Woodward Extracting Forceps	17
Scissors, Wire and Crown	100,101		
Seldin Root Elevators	28	Y	
Semken Tissue Forceps	102	Yankauer Suction Tubes	126
Senn-Miller Retractors	109	Yonger Good Universal Curettes	48
Sets, Dental Instruments	132,133	Young Pliers, Wire, Arch and Clasp Bending	89
Sliding Caliper	129		
Spatulas	128	Z	
Spatulas, Cement	43,44,128	Zahle Wax Knife	69
Spirit Lamp	122		
Spring Calipers	129		
Stadler Pliers	89		
Standard Extracting Forceps, American Pattern	16-24		
Standard Forceps, Dressing	102		
Standard Scalpel Handles	98		
Standard Scissors, Operating	99		
Sterilization	119-124		
Sterilizing Drums	123		
Sternberg Retractors, Cheek and Lip	113		
Stieglitz Splinter Forceps	106		
Stockes Rubber Dam Instruments	74		



SURGIFUSION
INTERNATIONAL

Contact: +92 313 4260504

Email: info@surgifusionintl.com

Website: www.surgifusionintl.com